

Video ToolBox™

Software for digital imaging

Display, capture, record, edit, label, measure, print, e-mail

Video Toolbox™ is the easiest way to display and capture video on your computer; desktop or laptop. Save a single image, a movie, a time lapse sequence, or capture the image directly into your favorite image-editing software for enhancement or annotation. Print or e-mail images without leaving the program.

Video Toolbox™ provides the tools you need to produce high quality images for live display or capture. Interactively adjust brightness, contrast and color settings. Zoom-in up to 4x without pixelation. Make precision measurements directly on the live image. Add labels, markers, even freehand drawings right on the live image.

Video Tools™ are plug-ins which add even more capability to Video ToolBox™. Plug-ins are available for image enhancement and automation. You can even write your own Video Tools.

Video ToolBox works with most USB, USB 2.0 and FireWire cameras which use DirectX WDM drivers. For analog video cameras, ask about our USB VideoLink adapters which put full-speed video into any computer including laptops.

Features:

- Display live video
- Capture images, time lapse sequences and movies.
- Capture images directly into your favorite image editor.
- Interactively adjust camera settings.
- Live zoom (without pixelation) up to 4x.
- Measurement on the live image using multiple measurement tools.
- Calibrate for different magnifications — store and recall up to 8 different calibrations.
- Add labels and markers over the live image. Drag to move.
- Place a time stamp anywhere on the image.
- Add an optional GPS receiver and your image will be marked with the exact location and time.
- Video Tool plug-ins provide additional capabilities for live image enhancement and automated measurement.

System Requirements:

Pentium II 500 MHz or faster (Pentium III+ recommended), Windows 2000/Me/XP, Internet Explorer 5.0 (or later) installed, USB Port, 128 MB of Free Memory, Video card with 8 MB of video RAM, minimum 800 x 600 resolution at 16-bit (65,000 colors); DirectDraw drivers recommended.

The image displays the Video Toolbox software interface with several windows overlaid on a live video feed of a cricket's back leg. The main window shows the video feed with a toolbar at the top and a status bar at the bottom. Overlaid windows include: 'Measure' showing a red circle and a red line with the word 'Sharp' written in red; 'Set Label' showing a text box with 'Cricket's back leg' and font settings; 'Zoom' showing a magnified view of the leg; and 'Video Averager' showing sliders for brightness, contrast, and gamma. The interface includes icons for Capture Mode, Contrast, Zoom, Measure, Draw, Label, Settings, Print, Help, Freeze, Live, and Capture.

Enhance with optional video tools

zarbeco, L.L.C.
creative imaging solutions
12 Ridge Road, Randolph, NJ 07869 USA
Phone: 973-989-2673 Fax: 973-989-2674
sales@zarbeco.com www.zarbeco.com